



HyperSpectral Imaging Microscope

User Guide





Table of Contents

| | |
|---|----|
| Chapter 1 Warning Symbol Definitions..... | 3 |
| Chapter 2 Safety..... | 4 |
| Chapter 3 Description..... | 5 |
| 3.1. Specifications Overview..... | 6 |
| Chapter 4 Getting Started..... | 7 |
| 4.1. General Setup..... | 7 |
| 4.2.1. Setting up the HYSPIM..... | 8 |
| 4.2.2. System Installation..... | 13 |
| Chapter 5 Capturing an Image..... | 16 |
| 5.1. Data Collection..... | 16 |



Chapter 1 Warning Symbol Definitions

Below is a list of warning symbols you may encounter in this manual or on your device.

| Symbol | Description |
|--------|---|
| | Direct Current |
| | Alternating Current |
| | Earth Ground Terminal |
| | On (Supply) |
| | Off (Supply) |
| | Caution: Risk of Electric Shock |
| | Caution: Burn hazard, Hot Surface, Do not touch |
| | Caution: Do Not Look Directly At Light |



Chapter 2 Safety

All statements regarding safety of operation and technical data in this user guide will only apply when the unit is operated correctly. Please read the following warnings and cautions carefully before operating the device



WARNING



Do not open the controller, electronic housing and detector module. There are no repairable parts in these products. Do not perform maintenance or service on the instrument unless specifically stated. Any modification or servicing of this system by unqualified personnel shall absolve of any liability. Only personnel authorized by and trained in the maintenance and repair of this equipment should remove any covers or enclosures or attempt any repairs or adjustments.



WARNING



The power supply supplied with this system is specifically designed for this system, do not replace or use in other devices.



Chapter 3 Description

The hyperspectral imaging microscope (HYSPIM) technique integrates the benefits of hyperspectral imaging (HSI) and microscopic imaging. HYSPIM offers both spatial and spectral information on substances, as well as their chemical composition at the molecular or cellular level. This makes HYSPIM a promising tool for nondestructive sample evaluation in various fields such as medicine, pathology, pharmaceuticals, life sciences, and the food industry.

The HYSPIM microscope system is designed with a special and precise optical system and an excellent mechanical system, which enables narrow spectral regions to be imaged individually.

A monochrome camera helps to obtain images with little spectral data for each pixel. Since the mechanical system has a very low error rate and does not move or swap during measurement, the data is not subject to pixel shifting or image registration issues.

The HYSPIM microscope comprises three key elements: a variable light source, a microscope, and software. The variable light source covers a spectral range from 400 nm to 1000 nm, with a spectral bandwidth of 6 nm. This source directs light to the microscope, positioned before the condenser lens. Here, the light undergoes reflection by a 45-degree mirror before illuminating the sample. Subsequently, a high-resolution camera integrated into the microscope captures the sample image.

The HYSPIM microscope is paired with sophisticated software that initiates the imaging process. This involves modifying the light source's wavelength while simultaneously recording the sample image, resulting in the construction of a data cube. The entire process seamlessly synchronizes with the software. Additionally, users can manually select a specific wavelength, facilitating the capture of a single image with a 6 nm bandwidth. This capability proves valuable for real-time examination of the sample under a specific narrow-band light source.

For users who prefer to utilize their existing microscope, the HYSPIM system provides the option to separately purchase the light source, software, and camera. For inquiries regarding this option, please reach out to our dedicated support department.





3.1. Specifications Overview

| SPECIFICATIONS | VALUE |
|-----------------------------|---|
| Spectral Range (nm) | Sensitive from 400 to 1000 nm |
| Spectral resolution | 6nm spectral resolution (100 Spectral channels) |
| Spatial resolution | limited by the microscope objective NA |
| Imaging method | Fast global mapping (snapshot) |
| Camera Type | CMOS Monochrome_Basler camera |
| Number of Active Pixels | 3800(H) x 2178 (V) |
| Pixel Size | 2.5 μm x 2.5 μm |
| Vertical Hardware Binning | Continuous Integer Values from 1 to 4 |
| Horizontal Hardware Binning | Continuous Integer Values from 1 to 4 |
| Exposure Time [us] | 51.0 to 1000 ms |
| Power Supply | Input : 230 V — 50 HZ Output : 24 V— 10 A |
| Halogen lamp | 24.0 V— 250 w |
| Color temperature | 3550 K |
| Weight | \approx 45 kg |
| Dimension | 85 cm \times 45 cm \times 45 cm (L,W,H) |
| Software | Python(windows 10&11) |



Chapter 4 Getting Started

4.1. General Setup



Controller Driver

12 V,1.5 A Power Supply

USB-3 Cable

Driver Cable

Power Cable

Monochrome Camera



4.2.1. Setting up the HYSPIM

1. Turn the microscope head to a suitable position



2. Remove the cap





3. Take the camera



4. Put the camera on the microscope and tighten the screw. In the next section there is an instruction to tune the direction of the camera





5. Plug the power cord to the light source power supply



6. Plug the cable from light source to the power supply



7. Plug the cable from driver box to the HYSPIM



8. Plug the cable from power supply adapter to the driver box





9. Connecting the communication cable of the HYSPIM controller to the driver box



10. Connecting the camera cable to the camera








4.2.2. System preparation

1. In the first step, you should adjust the central spot of light in the middle of the eyepiece of the microscope by means of the condenser adjustment screws





2. Then, using the height adjustment screw of the condenser, spread the light evenly so that the entire field of view of the objective is illuminated and the sample can be seen well.

**NOTE**

Note that all the above steps must be performed for different microscope objectives.



4. In the last step, you need to set the direction of the camera. First set the wavelength to 550 nm, then wait for the system to reach the selected wavelength. Second, using a 4x objective, try to see the line on the slide horizontally on the software screen







Chapter 5 Capturing an Image

5.1. Data Collection

1. Open the HSIM application on the computer
2. Select the Record Data from the Software window




3. Please select the driver control related to the device from the software's "com" button.
4. If you can't find the driver, click the Update button and try again. If the driver is still not found, check the device's connection control.






5. In the first step, you need to adjust the exposure time. First, set the wavelength to 550nm, then wait for the system to reach the selected wavelength. Second, make the necessary adjustments for exposure and gain




6. Now, by adjusting the Exposure level, it is possible to control the image intensity.






7. After adjusting the intensity, a) press the 'Return Home' button, b) change the position of the sample to a part where it is clear and not dirt exists, c) take the reference



WARNING


Note that before pressing the reference button, the exp value must match the value set at the time of adjusting the light with 550nm.

8. Press the Set Reference button.






9. Wait for the 'Reference' button on to turn from green to blue and for the checkmark to become active.




10. Press the 550nm button again and adjust the focus on the sample.






11. Now press the Return Home button.



12. Press the Start Scan button and wait for the data acquisition to complete.





13. After the window appears, you can save your data cube in the desired location.

