



Ratin Gene Poly A Polymerase Kit

Kit for addition of AMP from ATP to the 3' end of RNA

Components of the kit:

- RG Poly(A) Polymerase (CN: RG-Poly A), 50 μ L
- RG ATP 10mM (CN:RG-ATP), 100 μ L
- RG Poly(A) Polymerase Reaction Buffer 10X (CN: RG-Poly A Buffer), 100 μ L

*Keep at -20 °C

**Work on ice

Reaction Mix:

Poly(A) Polymerase Enzyme	0.5-1 μ L
Reaction Buffer 10X	2 μ L
ATP(10mM)	2 μ L
RNA Template	0.5-2 μ g
(RNase/DNase free Water)	variable
Total valume	20μL

- Incubate reaction at 37°C for 30 minutes.
- Stop the reaction by adding EDTA to a final concentration of 10 mM or directly proceeding to the cleanup step.



Ratin Gene Company

Keep Calm & Put Your Lab Coat On



0913 63 62 862

Ratingene.com

0913 63 62 963



Ratingenecompany@gmail.com