PrismaTech®

Benchtop Multimeter

Conductivity, Resistivity, Salinity, TDS

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

CONDUCTIVITY RANGE	0-500 ms/cm
RESISTIVITY RANGE	0-2000 Mohm*cm
SALINITY RANGE	0-70°
TDS RANGE	0-2000 mg/L
ACCURACY	± 0.5% Full Scale
TEMPERATURE COMPENSATION	Continues
ELECTRODS MATERIAL	AISI 316L / Gerafit / Platinum
ENCLOSURE MATERIAL	Anodized Aluminium
HEAD WETTED MATERIAL	PTFE
SAMPLE TEMPERATUE	-25-100°C
POWER SUPPLY	220VAC
PROTECTION	IP54
INTERNAL SENSORS	Board Voltages and Temperature
SERIAL OUTPUT	One Modbus RTU RS485
USB PORT	Availabe on HMI Display models
DATA LOGGER	10000 Samples (EC and Temperature) On HMI Display Model

Electrical Conductivity (EC), Electrical Resistivity, Salinity and TDS relate to the ability of a ionic solution to pass an electric current. In solutions the current is carried by cations and anions whereas in metals it is carried by electrons.

How well a solution conducts electricity depends on a number of factors

- Concentration
- Mobility of ions
- Valence of ions
- Temperature

In aqueous solutions the level of ionic strength varies from the low conductivity of ultra pure water to the high conductivity of concentrated chemical samples. The prob of the multimeter would be diped in the measuring sample and the parameters of the solution is calculated by applying an alternating electrical current (I) to two/four electrodes immersed in the solution and measuring the resulting voltage (V). During this process, the cations migrate to the negative electrode, the anions to the positive electrode and the solution acts as an electrical conductor.

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
PrismaTech[®] Benchtop Multimeter is powered by 220VAC without need for any external power supplier. There is also a USB port on the back of the device which the operator can attach a USB flash disk to it and save the stored data on it by pressing the "save as CSV" touch button and print the data directly on paper by attaching a USB cable to the port and the special printer.

The USB_B port which is placed beside the USB port on the back of the device can be used to connect the device to the computer and BPTM4_Lab Software.

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All setting and calibration operations of **PrismaTech®** Benchtop Multimeter can also be created through the BPTM4_Lab computer software. The computer can be easily connected to the sensor via the RS485 serial output and USB_Port so the device settings can be modified and the parameters can be monitored and saved to the computer from a distance up to 1000meters far from the device.

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- ✓ Abrupt measurement and high accuracy
- ✓ Including RS485 port with Modbus protocol to communicate with a computer up to a 1000m far
- ✓ Creating alarm in bad working conditions or error in communication, decreasing voltage etc.
- ✓ Isolated against electrical noise and mechanical vibrations
- ✓ Applicable to measure the conductivity of a wide range of solutions
- ✓ With high resolution touch HMI 4.3" display
- ✓ Ability to store measured data on the device memory and save them on a USB flash disk
- ✓ Ability to connect directly to the printer and print the results on the paper
- ✓ 220VAC Power supply without need to adaptor
- ✓ User friendly software and easy setting and calibration procedure
- ✓ One year warranty and ten-year real after-sales service

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✓ Measurement of the four parameters^{*} of different electrolyte solutions.

✓ Measurement of the four parameters of solutions in every station of the purification process

✓ Measurement of the four parameters of acids, minerals, ammonia, organic solutions, salty solutions etc.

✓ Measurement of the four parameters of many chemical solutions derived from oil, slurry, Hydrogen peroxide, hydrofluoric acids etc

✓ Measurement of the concentration of vinegar and salt water in pickles and dairy industry

✓ Measurement of the concentration of hydrofluoric acid in acid washing lines in the steel rolling factories

✓ Measurement of the four parameters and water hardness in the cooling systems and boilers in the power plants

✓ Measurement of the four parameters of solutions of the edible and industrial salts

...